


SU ARVID

Multitalent in usage for Grain, Feeding and Biogas



Advantages:

- very high and stable yield due to a longer grain filling phase
- 2018 drought: European winner in yield
- short growing type with good lodging tolerance and still high dry matter
- one of the best double use variety in SU Portfolio (whole crop and grain)
- Above average in leaf health

Cultivation:

Highest flexibility for your use of rye and no limitation for certain soil types. SU Arvid shows mid-high enzyme activity which is needed in years with very high falling numbers. In difficult years SU Arvid should be harvested early to ensure stable falling numbers.

Short profile:

---- = very low resistance/early/short,

++++ = very high resistance/late/long

Plant Development

Maturity




Vitality

Lodging resistance



Health

Ergot



Quality

Hagberg / Falling number



Raw protein content



SU ARVID

Multitalent in usage for Grain, Feeding and Biogas

Plant development and vitality:

Ear emergence	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Maturity	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>

Yield structure:

Ears/m ²	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Grains/ear	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
TGW	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>

Straw:

Plant height	short/medium plant height								
Lodging resistance	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>

Disease resistance:

Mildew	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Rhynchosporium	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Brown rust	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Ergot	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>

Quality:

Hagberg / Falling number	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Raw protein content	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>
Amylogram-values	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>	<div></div>